

ZytoLight® SPEC YWHAE Dual Color Break Apart Probe




Background

The ZytoLight® SPEC YWHAE Dual Color Break Apart Probe (PL134) is intended to be used for the qualitative detection of translocations involving the human YWHAE gene at 17p13.3 in formalin-fixed, paraffin-embedded specimens, such as endometrial stromal sarcoma (ESS), by fluorescence *in situ* hybridization (FISH). The probe is intended to be used in combination with the ZytoLight® FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20). The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel. The probe is intended to be used as an aid to the differential diagnosis of ESS and therapeutic measures should not be initiated based on the test result alone.


Probe Description

The ZytoLight® SPEC YWHAE Dual Color Break Apart Probe is composed of:

- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~10 ng/µl), which target sequences mapping in 17p13.3** (chr17:1,339,752-1,716,668) proximal to the YWHAE breakpoint region.
- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/µl), which target sequences mapping in 17p13.3** (chr17:791,171-1,124,746) distal to the YWHAE breakpoint region.
- Formamide based hybridization buffer




Ideogram of chromosome 17 indicating the hybridization locations.



SPEC YWHAE Probe map (not to scale).

Results

In an interphase nucleus of a normal cell lacking a translocation involving the 17p13.3 band, two orange/green fusion signals are expected representing two normal (non-rearranged) 17p13.3 loci. A signal pattern consisting of one orange/green fusion signal, one orange signal, and a separate green signal indicates one normal 17p13.3 locus and one 17p13.3 locus affected by a translocation.



Endometrial stromal sarcoma tissue section with translocation affecting the YWHAE gene as indicated by one non-rearranged orange/green fusion signal, one orange, and one separate green signal.

Prod. No.	Product	Label	Tests* (Volume)
Z-2175-50	ZytoLight SPEC YWHAE Dual Color Break Apart Probe CE 0124 IVD	●/●	5 (50 µl)
Related Products			
Z-2028-5	ZytoLight FISH-Tissue Implementation Kit CE IVD		5
Incl. Heat Pretreatment Solution Citric, 150 ml; Pepsin Solution, 1 ml; Wash Buffer SSC, 210 ml; 25x Wash Buffer A, 50 ml; DAPI/DuraTect-Solution, 0.2 ml			

* Using 10 µl probe solution per test. IVD labeled products are only available in certain countries. All other countries research use only! Please contact your local dealer for more information.

**According to Human Genome Assembly GRCh37/hg19