



## FlexISH- Tissue Implementation Kit

**REF** Z-2182-5 Σ 5

**REF** Z-2182-20 Σ 20

For use in fluorescence *in situ* hybridization (FISH)  
procedures

4250380N8486



In vitro diagnostic medical device  
according to IVDR (EU) 2017/746

### 1. Intended use

The FlexISH-Tissue Implementation Kit is intended to be used in combination with FlexISH probes on formalin-fixed, paraffin-embedded specimens by fluorescence *in situ* hybridization (FISH).

The product is intended for professional use only. All tests using the product should be performed in a certified, licensed anatomic pathology laboratory under the supervision of a pathologist/human geneticist by qualified personnel.

### 2. Test principle

The fluorescence *in situ* hybridization (FISH) technique allows for the detection and visualization of specific nucleic acid sequences in cell preparations. Fluorescently-labeled DNA fragments, so called FISH probes, and their complementary target DNA strands in the preparations are co-denatured and subsequently allowed to anneal during hybridization. Afterwards, unspecific and unbound probe fragments are removed by stringency washing steps. After counterstaining the DNA with DAPI, hybridized probe fragments are visualized using a fluorescence microscope equipped with excitation and emission filters specific for the fluorochromes with which the FISH probe fragments have been directly labeled.

### 3. Reagents provided

The FlexISH-Tissue Implementation Kit is available in two sizes and is composed of:

Code	Component	Quantity		Container
		20	5	
PT1	<u>Heat Pretreatment Solution Citric</u>	500 ml	150 ml	Screw-cap bottle (large)
ES1	<u>Pepsin Solution</u>	4 ml	1 ml	Dropper bottle, white cap
WB10	<u>5x FlexISH Wash Buffer</u>	500 ml	150 ml	Screw-cap bottle (large)
MT7	<u>DAPI/DuraTect-Solution</u>	0.8	0.2	Reaction vessel, blue lid
	Instructions for use	1	1	

**Z-2182-5 (5 tests):** Components **ES1** and **MT7** are sufficient for 5 reactions. Component **WB10** is sufficient for 3x 3 staining jars of 70 ml each. Component **PT1** is sufficient for 2 staining jars of 70 ml each.

**Z-2182-20 (20 tests):** Components **ES1** and **MT7** are sufficient for 20 reactions. Component **WB10** is sufficient for 1x 3 staining jars of 70 ml each. Component **PT1** is sufficient for 7 staining jars of 70 ml each.

### 4. Materials required but not provided

- FlexISH probe
- Positive and negative control specimens
- Microscope slides, positively charged
- Water bath (37°C, 98°C)
- Hot plate or hybridizer
- Humidity chamber + hybridization oven or hybridizer
- Adjustable pipettes (10 µl, 30 µl)
- Staining jars or baths
- Timer
- Calibrated thermometer
- Ethanol or reagent alcohol
- Xylene
- Deionized or distilled water
- Coverslips (22 mm x 22 mm, 24 mm x 60 mm)
- Rubber cement, e.g., Fixogum Rubber Cement (Prod. No. E-4005-50/-125) or similar
- Adequately maintained fluorescence microscope (400-1000x)
- Immersion oil approved for fluorescence microscopy
- Appropriate filter sets

### 5. Storage and handling

Store at 2-8 °C in an upright position. Additionally, the DAPI/DuraTect-Solution (**MT7**) must be stored protected from light. Return to storage conditions immediately after use. Do not use reagents beyond expiry date indicated on the label. The product is stable until expiry date indicated on the label when handled accordingly.

### 6. Warnings and precautions

- Read the instructions for use prior to use!
- Do not use the reagents after the expiry date has been reached!
- This product contains substances (in low concentrations and volumes) that are harmful to health. Avoid any direct contact with the reagents. Take appropriate protective measures (use disposable gloves, protective glasses, and lab garments)!
- Report any serious incident that has occurred in relation to the product to the manufacturer and the competent authority according to local regulations!

- If reagents come into contact with skin, rinse skin immediately with copious amounts of water!
- A material safety data sheet is available on request for the professional user.
- Do not reuse reagents, unless reuse is explicitly permitted!
- Avoid cross-contamination of samples as this may lead to erroneous results.
- The specimens must not be allowed to dry during the hybridization and washing steps.
- DAPI/DuraTect-Solution (MT7) should not be exposed to light, especially strong light, for a longer period of time, i.e., all steps should be accomplished, where possible, in the dark and/or using lightproof containers!

#### Special labelling of ES1:

EUH208	Contains Pepsin A. May produce an allergic reaction.
EUH210	Safety data sheet available on request.

#### Hazard and precautionary statements PT1 and WB10:

The hazard determining component is a reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6] (3:1).



#### Warning

H317	May cause an allergic skin reaction.
P261	Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves/ protective clothing/ eye protection/ face protection.
P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice/attention.
P362+P364	Take off contaminated clothing and wash it before reuse.

#### Hazards and precautionary statements for MT7:

This product is not classified as hazardous according to Regulation (EC) No. 1272/2008.

#### 7. Limitations

- For *in vitro* diagnostic use.
- For professional use only.
- For non-automated use only.
- The clinical interpretation of any positive staining, or its absence, must be done within the context of clinical history, morphology, other histopathological criteria as well as other diagnostic tests. It is the responsibility of a qualified pathologist/human geneticist to be familiar with the ISH probes, reagents, diagnostic panels, and methods used to produce the stained preparation. Staining must be performed in a certified, licensed laboratory under the supervision of a pathologist/human geneticist who is responsible for reviewing the stained slides and assuring the adequacy of positive and negative controls.
- Specimen staining, especially signal intensity and background staining, is dependent on the handling and processing of the specimen prior to staining. Improper fixation, freezing, thawing, washing, drying, heating, sectioning, or contamination with other specimens or fluids may produce artefacts or false results. Inconsistent results may result from variations in fixation and embedding methods, as well as from inherent irregularities within the specimen.

- The performance was validated using the procedures described in the instruction for use of the respective ZytoVision probe and implementation kit. Modifications to these procedures might alter the performance and have to be validated by the user. This IVD is only certified as CE when used as described in this instruction for use within the scope of the intended use.

#### 8. Interfering substances

Blood cells present in the specimen might exhibit autofluorescence which hinders signal recognition.

The following fixatives are incompatible with ISH:

- Bouin's fixative
- B5 fixative
- Acidic fixatives (e.g., picric acid)
- Zenker's fixative
- Alcohols (when used alone)
- Mercuric chloride
- Formaldehyde/zinc fixative
- Hollande's fixative
- Non-buffered formalin

#### 9. Preparation of specimens

Recommendations:

- Fixation in 10% neutrally buffered formalin for 24 h at room temperature (18°C-25°C).
- Sample size  $\leq 0.5 \text{ cm}^3$ .
- Use premium quality paraffin.
- Embedding should be carried out at temperatures lower than 65°C.
- Prepare 2-4  $\mu\text{m}$  microtome sections.
- Use positively charged microscope slides.
- Fix for 2-16 h at 50-60°C.

#### 10. Preparatory treatment of the device

5x FlexISH Wash Buffer (WB10) is to be pretreated according to the instructions in 11. "Assay procedure". All other kit reagents are ready-to-use. No reconstitution, mixing, or dilution is required.

#### 11. Assay procedure

##### 11.1 Day 1

##### Preparatory steps

1. *Preparation of two ethanol series (70%, 90%, and 100% ethanol solutions):* Dilute 100% ethanol with deionized or distilled water. These solutions can be stored in suitable containers and can be re-used for up to 160 slides.
2. *Heat Pretreatment Solution Citric (PT1):* Fill a staining jar and warm to 98°C.
3. *FlexISH Probe:* Bring to room temperature before use, protect from light. Prior to opening the vial, mix by vortexing and spin down briefly.

##### Pretreatment (dewax/proteolysis)

1. Incubate slides for 2x 5 min in xylene.
2. Incubate in 100%, 100%, 90%, and 70% ethanol, each for 2 min.
3. Wash 2x 2 min in deionized or distilled water.
4. Incubate for 20 min in pre-warmed Heat Pretreatment Solution Citric (PT1) at 98°C.

*We recommend not to use more than eight slides per staining jar. After immersing the slides, check the temperature of the Heat Pretreatment Solution Citric inside the jar and start time as soon as the temperature of the solution has reached at least 95°C.*

5. Transfer slides immediately to deionized or distilled water, wash for 2x 2 min and drain off or blot off the water.
6. Apply (dropwise) Pepsin Solution (ES1) to the specimens and incubate for 15 min at 37°C in a humidity chamber.

**ES1** may form precipitates, which do not affect the quality.

Depending on multiple factors, e.g., nature and duration of fixing, thickness of sections, and nature of specimens, different incubation times may be required. As an incubation guideline, we recommend an incubation time of 2-30 min. As a general rule, we recommend to ascertain the optimum time for proteolysis in pre-tests.

7. Wash for 2x 2 min in deionized or distilled water.
8. Dehydration: in 70%, 90%, and 100% ethanol, each for 1 min.
9. Air dry sections.

Make sure to completely dry sections prior to probe application since residual moisture may reduce signal intensity and/or affect specimen morphology.

**Denaturation and hybridization**

1. Pipette 10 µl of the FlexISH Probe onto each pretreated specimen.  
*Avoid long exposure of the probe to light.*
2. Cover specimens with a 22 mm x 22 mm coverslip (avoid trapped bubbles) and seal the coverslip.

*We recommend using rubber cement (e.g., Fixogum) for sealing.*

3. Place slides on a hot plate or hybridizer and denature specimens for 10 min at 75°C.
4. Perform hybridization for 2 h up to 16 h (i.e. overnight) at 37°C by either transferring the slides to a hybridizer or to a humidity chamber and a hybridization oven.

*It is essential that specimens do not dry out during the hybridization step.*

**11.2 Day 1 or Day 2**

**Preparatory steps**

1. Preparation of 1x FlexISH Wash Buffer: Dilute 1 part 5x FlexISH Wash Buffer (WB10) with 4 parts deionized or distilled water. Fill three staining jars with the 1x FlexISH Wash Buffer, pre-warm one jar to 72°C and keep two jars at room temperature.
2. DAPI/DuraTect-Solution (MT7): Bring to room temperature before use, protect from light.

**Post-hybridization and detection**

1. Carefully remove the rubber cement or glue.
2. Remove the coverslips by submerging in 1x FlexISH Wash Buffer at room temperature for 1-2 min.

*To facilitate the removal of the coverslip, this step can alternatively be performed for 2 min at 37°C.*

3. Wash using 1x FlexISH Wash Buffer for 10 min at 72°C.  
*The 1x FlexISH Wash Buffer should be pre-warmed. Check with a thermometer if necessary. Do not use more than eight slides per staining jar.*
4. Wash using 1x FlexISH Wash Buffer for 3 min at room temperature.
5. Incubate the slides in 70%, 90%, and 100% ethanol, each for 1 min.
6. Air dry the samples protected from light.
7. Pipette 25 µl DAPI/DuraTect-Solution (MT7) onto the slides. Avoiding trapped bubbles, cover the samples with a coverslip (24 mm x 60 mm). Incubate in the dark for 15 min.

*Using a pipette tip which has been cut off to increase the size of the opening, can make the pipetting process easier. Avoid long exposure to light.*

8. Store the slides in the dark. For longer storage periods, this should take place at 2-8°C.
9. Evaluation of the sample material is carried out by fluorescence microscopy. Filter sets for the following wavelength ranges are required:

Fluorescent dye	Excitation	Emission
ZyBlue	418 nm	467 nm
ZyGreen	503 nm	528 nm
ZyGreen 2.0	493 nm	518 nm
ZyGold	532 nm	553 nm
ZyOrange	547 nm	572 nm
ZyRed	580 nm	599 nm

**12. Interpretation of results**

With the use of appropriate filter sets in interphases or metaphases of normal cells or cells without aberrations of chromosomes, two signals per probe/fluorescence label appear, except for probes targeting X and/or Y chromosomes, resulting in none to two signals per probe/fluorescence label, depending on the gender. In cells with chromosomal aberrations, a different signal pattern can be visible in interphases or metaphases. For more details on the interpretation of results, please refer to the respective probe manual.

**13. Recommended quality control procedures**

Refer to the instructions for use of the respective ZytoVision probe.

**14. Performance characteristics**

Refer to the instructions for use of the respective ZytoVision probe.

**15. Disposal**

The disposal of reagents must be carried out in accordance with local regulations.

**16. Troubleshooting**

Any deviation from the operating instructions can lead to inferior staining results or to no staining at all. Please refer to [www.zytovision.com](http://www.zytovision.com) for more information.

**Weak signals or no signals at all**

Possible cause	Action
Specimen has not been properly fixed	Optimize fixing time and fixative
Proteolytic pretreatment not carried out properly	Optimize pepsin incubation time, increase or decrease if necessary
Probe evaporation	When using a hybridizer, the use of the wet stripes/water filled tanks is mandatory. When using a hybridization oven, the use of a humidity chamber is required. In addition, the coverslip should be sealed completely, e.g., with Fixogum, to prevent drying-out of the sample during hybridization.
Inappropriate filter sets used	Use filter sets appropriate for the fluorochromes of the probe. <i>Triple-bandpass filter sets provide less light compared to single or dual-bandpass filter sets. Consequently, the signals may appear fainter using these triple-bandpass filter sets.</i>

**Cross hybridization signals; noisy background**

Possible cause	Action
Incomplete dewaxing	Use fresh solutions; check duration of dewaxing
Proteolytic pretreatment too strong	Reduce pepsin incubation time
Slides cooled to room temperature before hybridization	Transfer the slides quickly to 37 °C

**Morphology degraded**

Possible cause	Action
Specimen has not been properly fixed	Optimize fixing time and fixative
Proteolytic pretreatment not carried out properly	Optimize pepsin incubation time, decrease if necessary
Insufficient drying before probe application	Extend air-drying

**Overlapping nuclei**

Possible cause	Action
Inappropriate thickness of tissue sections	Prepare 2-4 $\mu\text{m}$ microtome sections

**Specimen floats off the slide**

Possible cause	Action
Proteolytic pretreatment too strong	Reduce pepsin incubation time

**Weak counterstain**

Possible cause	Action
Low concentrated DAPI solution	Use <u>DAPI/DuraTect-Solution (ultra)</u> (Prod. No. MT-0008-0.8) instead
DAPI incubation time too short	Adjust DAPI incubation time

**17. Literature**

- Kievits T, et al. (1990) *Cytogenet Cell Genet* 53: 134-6.
- Wilkinson DG: In Situ Hybridization, A Practical Approach, *Oxford University Press* (1992) ISBN 0 19 963327 4.

**18. Revision**

Revision	Description of the change
1.2.1	11. Assay procedure ZyGreen 2.0 added



[www.zytovision.com](http://www.zytovision.com)

Please refer to [www.zytovision.com](http://www.zytovision.com) for the most recent instructions for use as well as for instructions for use in different languages.

Our experts are available to answer your questions.  
Please contact [helptech@zytovision.com](mailto:helptech@zytovision.com)



ZytoVision GmbH  
Fischkai 1  
27572 Bremerhaven/ Germany  
Phone: +49 471 4832-300  
Fax: +49 471 4832-509  
[www.zytovision.com](http://www.zytovision.com)  
Email: [info@zytovision.com](mailto:info@zytovision.com)

**Trademarks:**  
ZytoVision® and F/exSH® are trademarks of ZytoVision GmbH.