

VisionArray® Analysis Package

For visualization and quick analysis of the VisionArray® Chips data

Legend

- 1 Scanner 8100
- 2 Hand scanner
- 3 Slide Holder
- 4 Laptop
- 5 USB-Hub
- 6 External hard drive
- 7 Computer mouse



Prod. No. Product

E-4060-1

VisionArray Analysis Package CE [IVD]

Incl. Scanner 8100; Slide Holder; Hand Scanner; PC with preinstalled VisionArray Analyzer Software; USB-Hub; External Hard Drive; Computer Mouse

VisionArray® Analyzer Software

- Simple visualization and quick analysis of the VisionArray® Chip data
- Analysis of a chip and the report of the results can be achieved in just a few minutes
- Program navigation is easy and intuitive for the user
- Scans are stored including all sample and chip data in an integrated database on the enclosed external hard drive

VisionArray® Analyzer
28/12/2016 10:49:37

Start Archive x 5084_A001-PK06-094 x

Scan-ID: 5084	Name:	Slide-ID: A001-PK06-094	Date of Scan: 28/12/2016 10:46
Case No.: HPV Testing	Date of Birth:	Chip-Lot: PK06	Chip File: VisionArray HPV Chip 1.0

Zoom: 13% | Overlay On/Off

Expert Mode:

✓ Hybridization/Grid ✓ Positive Control
✓ Duplicates

Target	Intensity 1	Intensity 2
HPV11	50.08	29.71%
HPV18	139.04	82.50%
HPV82 (IS39)	70.16	41.62%
Hyb./Grid	168.55	100.00%
Pos.	173.20	102.76%

Time Stamp	Activity	De
28/12/2016 10:47:53	Scan created	
28/12/2016 10:47:53	Automatic Analysis	

Next

CE [IVD] only available in certain countries. All other countries research use only! Please contact your local dealer for more information.